




**BAKERY**



**81601**

**1 Cell Tray**

**O.D.:** 2.750 x 2.750 x 0.250

**I.D.:** 2.500 x 2.500 x 0.250

**Mat.:** OPS, HIPS



**2656**

**1 Cell Tray**

**O.D.:** 3.125 x 2.625 x 1.750

**I.D.:** 2.375 x 1.500 x 1.750

**Mat.:** OPS, HIPS



**2103**

**1 Cell Tray**

**O.D.:** 3.250 x 3.000 x 1.375

**I.D.:** 2.750 x 2.000 x 1.375

**Mat.:** OPS, HIPS



**86375**

**1 Cell Tray**

**O.D.:** 3.500 x 2.375 x 0.750

**I.D.:** 3.000 x 1.875 x 0.625

**Mat.:** OPS, HIPS



**86322**

**1 Cell Tray**

**O.D.:** 3.750 x 2.750 x 2.000

**I.D.:** 3.250 x 2.250 x 2.000

**Mat.:** OPS, HIPS



**2414**

**1 Cell Tray**

**O.D.:** 4.125 x 2.750 x 0.875

**I.D.:** 3.875 x 2.250 x 0.875

**Mat.:** OPS, HIPS



**3501**

**1 Cell Tray**

**O.D.:** 4.125 x 2.750 x 1.500

**I.D.:** 3.375 x 2.250 x 1.500

**Mat.:** OPS, HIPS



**2415**

**1 Cell Tray**

**O.D.:** 4.750 x 2.750 x 0.875

**I.D.:** 4.125 x 2.250 x 0.875

**Mat.:** OPS, HIPS



**82980**

**1 Cell Tray**

**O.D.:** 4.750 x 4.750 x 1.000

**I.D.:** 4.500 x 4.500 x 1.000

**Mat.:** OPS, HIPS



**3007**

**Pie Wedge, Deep**

**O.D.:** 5.000 x 4.000 x 1.750

**I.D.:** 4.500 x 3.500 x 1.500

**Mat.:** OPS, HIPS



**3013**

**Pie Wedge, Shallow**

**O.D.:** 5.000 x 4.000 x 1.250

**I.D.:** 3.375 x 2.875 x 1.325

**Mat.:** OPS, HIPS



**88401**

**1 Cell Tray**

**O.D.:** 5.000 x 4.000 x 1.625

**I.D.:** 4.500 x 3.500 x 1.625

**Mat.:** OPS, HIPS



**85353**

**1 Cell Tray**

**O.D.:** 5.250 x 3.125 x 1.500

**I.D.:** 4.875 x 2.250 x 1.500

**Mat.:** OPS, HIPS



**87423**

**1 Cell Tray**

**O.D.:** 5.375 x 5.125 x 1.125

**I.D.:** 4.125 x 4.125 x 1.000

**Mat.:** OPS, HIPS




**81672**

**1 Cell Tray**

**O.D.:** 6.125 x 4.500 x 2.125

**I.D.:** 5.500 x 3.750 x 2.125

**Mat.:** OPS, HIPS



**85357**

**1 Cell Tray**

**O.D.:** 6.375 x 4.750 x 1.625

**I.D.:** 5.875 x 4.125 x 1.625

**Mat.:** PETE



**2104**

**1 Cell Tray**

**O.D.:** 6.500 x 3.000 x 1.500

**I.D.:** 6.000 x 1.750 x 1.500

**Mat.:** PETE




**87301**

**1 Cell Tray**

**O.D.:** 6.625 x 5.875 x 1.375

**I.D.:** 6.000 x 5.250 x 1.375

**Mat.:** OPS, HIPS



**2201**

**1 Cell Tray**

**O.D.:** 6.750 x 4.750 x 1.000

**I.D.:** 6.000 x 4.125 x 0.875

**Mat.:** OPS, HIPS



**86101**

**1 Cell Tray**

**O.D.:** 6.750 x 2.625 x 1.125

**I.D.:** 6.250 x 2.250 x 1.125

**Mat.:** OPS, HIPS




**2606**

**1 Cell Tray**

**O.D.:** 7.000 x 7.000 x 1.500

**I.D.:** 6.500 x 6.500 x 1.500

**Mat.:** OPS, HIPS




**86345**

**1 Cell Tray**

**O.D.:** 7.250 x 2.750 x 2.000

**I.D.:** 7.625 x 2.250 x 2.000

**Mat.:** OPS, HIPS



**5701**

**1 Cell Tray**

**O.D.:** 7.375 x 5.250 x 0.875

**I.D.:** 7.000 x 4.750 x 0.875

**Mat.:** OPS, HIPS




**3015**

**1 Cell Tray**

**O.D.:** 7.625 x 5.250 x 0.500

**I.D.:** 6.750 x 4.500 x 0.375

**Mat.:** OPS, HIPS




**86012**

**1 Cell Tray**

**O.D.:** 7.625 x 2.375 x 2.000

**I.D.:** 7.500 x 1.875 x 2.000

**Mat.:** OPS, HIPS



**2352**

**1 Cell Tray**

**O.D.:** 8.000 x 4.375 x 2.000

**I.D.:** 7.875 x 3.625 x 1.875

**Mat.:** OPS, HIPS



**3525**


**1 Cell Tray**

**O.D.:** 8.000 x 4.000 x 1.625

**I.D.:** 7.500 x 3.125 x 1.625

**Mat.:** OPS, HIPS






**86501**

**1 Cell Tray**

**O.D.:** 8.125 x 5.000 x 1.000

**I.D.:** 7.875 x 4.250 x 1.000

**Mat.:** OPS, HIPS




**84501**

**1 Cell Tray**

**O.D.:** 8.250 x 3.250 x 0.875

**I.D.:** 8.000 x 2.750 x 0.875

**Mat.:** OPS, HIPS




**2200**

**1 Cell Tray**

**O.D.:** 8.375 x 4.000 x 1.375

**I.D.:** 7.875 x 3.500 x 1.250

**Mat.:** OPS, HIPS



**2106**

**1 Cell Tray**

**O.D.:** 8.625 x 7.250 x 0.675

**I.D.:** 8.250 x 6.750 x 0.500

**Mat.:** OPS, HIPS



**89801**

**Pie Liner (Fits 8 1/2" Pie Shell)**

**O.D.:** 8.750 x 8.750 x 1.000

**I.D.:** 6.500 x 6.500 x 1.000

**Mat.:** OPS, HIPS




**86391**

**1 Cell Tray, Sealable**

**O.D.:** 9.000 x 7.000 x 2.500

**I.D.:** 8.250 x 6.000 x 2.250

**Mat.:** PETE



**2220**

**1 Cell Tray**

**O.D.:** 9.250 x 5.375 x 1.500

**I.D.:** 8.500 x 4.500 x 1.500

**Mat.:** OPS, HIPS




**2113**

**1 Cell Tray**

**O.D.:** 9.500 x 4.000 x 0.625

**I.D.:** 9.200 x 3.500 x 0.675

**Mat.:** OPS, HIPS



**2108**

**1 Cell Tray**

**O.D.:** 9.625 x 7.250 x 0.625

**I.D.:** 9.250 x 7.375 x 0.500

**Mat.:** OPS, HIPS



**2221**

**1 Cell Tray**

**O.D.:** 9.625 x 4.875 x 1.750

**I.D.:** 9.375 x 4.125 x 1.375

**Mat.:** OPS, HIPS



**85351**

**1 Cell Tray**

**O.D.:** 9.625 x 6.500 x 1.625

**I.D.:** 9.250 x 5.750 x 1.625

**Mat.:** OPS, HIPS



**2351**

**1 Cell Tray**

**O.D.:** 10.000 x 5.500 x 1.375

**I.D.:** 9.125 x 4.875 x 1.375

**Mat.:** OPS, HIPS



**5001**

**1 Cell Tray**

**O.D.:** 10.375 x 3.750 x 1.500

**I.D.:** 10.000 x 3.000 x 1.500

**Mat.:** OPS, HIPS




**80202**

**1 Cell Tray, Sealable**

**O.D.:** 10.375 x 8.250 x 1.250

**I.D.:** 9.125 x 7.000 x 1.250

**Mat.:** OPS, HIPS, PETE




**2065**

**1 Cell Tray**

**O.D.:** 10.500 x 3.875 x 1.675

**I.D.:** 10.250 x 3.250 x 1.500

**Mat.:** OPS, HIPS




**2800**

**1 Cell Tray**

**O.D.:** 10.750 x 7.500 x 1.000

**I.D.:** 10.000 x 6.875 x 1.000

**Mat.:** OPS, HIPS



**5700**

**1 Cell Tray**

**O.D.:** 10.750 x 5.500 x 1.250

**I.D.:** 10.500 x 5.125 x 1.250

**Mat.:** OPS, HIPS



**3499**

**1 Cell Tray**

**O.D.:** 11.000 x 5.000 x 1.375

**I.D.:** 10.375 x 4.625 x 1.375

**Mat.:** OPS, HIPS



2025

2 Cell Tray

O.D.: 5.125 x 3.875 x 1.750

I.D.: 2.500 x 2.000 x 1.500

Mat.: OPS, HIPS



86802

2 Cell Tray

O.D.: 5.250 x 4.500 x 2.250

I.D.: 4.000 x 2.125 x 2.250

Mat.: OPS, HIPS



87231

2 Cell Tray

O.D.: 5.500 x 4.750 x 1.500

I.D.: 5.000 x 1.750 x 1.500

Mat.: PETE



2657

2 Cell Tray

O.D.: 6.750 x 5.125 x 1.125

I.D.: 6.250 x 2.250 x 1.125

Mat.: OPS, HIPS



2668

2 Cell Rectangle Tray

O.D.: 7.125 x 3.625 x 2.000

I.D.: 3.000 x 3.000 x 1.875

Mat.: OPS, HIPS



82775

2 Cell Tray

O.D.: 7.125 x 3.625 x 2.500

I.D.: 2.750 x 2.750 x 2.500

Mat.: PETE



2422

2 Cell Tray

O.D.: 10.000 x 5.000 x 1.250

I.D.: 4.500 x 4.500 x 1.250

Mat.: OPS, HIPS



80050

2 Cell Tray

O.D.: 10.000 x 5.500 x 2.250

I.D.: 4.375 x 5.000 x 2.250

4.375 x 5.000

Mat.: OPS, HIPS



87422

2 Cell Tray

O.D.: 10.250 x 5.250 x 1.000

I.D.: 4.000 x 4.000 x 1.000

Mat.: PETE



**80503**

**3 Cell Tray**

**O.D.:** 6.375 x 5.375 x 2.250

**I.D.:** 1.750 x 4.875 x 2.250

**Mat.:** OPS, HIPS



**2051**

**3 Cell Tray**

**O.D.:** 6.675 x 3.875 x 2.125

**I.D.:** 6.250 x 3.500 x 2.125

**Mat.:** OPS, HIPS



**2053**

**3 Cell Tray**

**O.D.:** 6.750 x 4.375 x 1.750

**I.D.:** 6.375 x 3.750 x 1.675

**Mat.:** OPS, HIPS



**2601**

**3 Cell Tray**

**O.D.:** 7.750 x 4.000 x 2.250

**I.D.:** 2.000 x 3.500 x 2.250

**Mat.:** OPS, HIPS



**85030**

**3 Cell Tray**

**O.D.:** 7.750 x 6.250 x 1.750

**I.D.:** 7.750 x 1.750 x 1.750

**Mat.:** OPS, HIPS



**2357**

**3 Cell Tray**

**O.D.:** 7.875 x 3.500 x 2.000

**I.D.:** 1.875 x 2.750 x 2.000

**Mat.:** OPS, HIPS



**2050**

**3 Cell Tray**

**O.D.:** 8.250 x 3.875 x 1.750

**I.D.:** 2.500 x 3.250 x 1.750

**Mat.:** OPS, HIPS



**2409**

**3 Cell Tray**

**O.D.:** 8.750 x 3.750 x 2.000

**I.D.:** 2.375 x 2.875 x 1.875

**Mat.:** OPS, HIPS



**88200**

**3 Cell Tray**

**O.D.:** 8.875 x 3.875 x 2.000

**I.D.:** 2.500 x 3.375 x 1.875

**Mat.:** OPS, HIPS





2190

3 Cell Tray

O.D.: 9.000 x 3.250 x 1.375

I.D.: 1.750 x 2.750 x 1.250

Mat.: OPS, HIPS



2195

3 Cell Tray

O.D.: 9.000 x 3.250 x 2.250

I.D.: 1.750 x 2.750 x 2.250

Mat.: OPS, HIPS



2403

3 Cell Tray

O.D.: 9.000 x 3.500 x 1.375

I.D.: 1.500 x 3.125 x 1.375

Mat.: OPS, HIPS



2410

3 Cell Tray

O.D.: 9.000 x 3.250 x 2.250

I.D.: 2.675 x 2.675 x 2.125

Mat.: OPS, HIPS



86203

3 Cell Tray

O.D.: 9.375 x 3.375 x 1.250

I.D.: 3.250 x 2.500 x 1.250

Mat.: OPS, HIPS



85017

3 Cell Tray

O.D.: 9.375 x 6.000 x 1.000

I.D.: 8.625 x 1.875 x 1.000

Mat.: OPS, HIPS



2411

3 Cell Tray

O.D.: 9.625 x 3.625 x 1.750

I.D.: 3.000 x 3.000 x 1.750

Mat.: OPS, HIPS



2434

3 Cell Tray

O.D.: 9.750 x 3.500 x 1.500

I.D.: 3.000 x 3.000 x 1.375

Mat.: OPS, HIPS



81344

3 Cell Tray

O.D.: 9.820 x 3.525 x 1.250

I.D.: 3.000 x 3.000 x 1.250

Mat.: OPS, HIPS



2412

3 Cell Tray

O.D.: 10.750 x 3.750 x 1.625

I.D.: 3.125 x 3.125 x 1.625

Mat.: OPS, HIPS



2413

3 Cell Tray

O.D.: 10.750 x 3.750 x 1.625

I.D.: 3.250 x 3.250 x 1.625

Mat.: OPS, HIPS



85040

4 Cell Tray

O.D.: 8.000 x 3.500 x 1.500

I.D.: 1.625 x 3.000 x 1.500

Mat.: OPS, HIPS



2424

4 Cell Tray

O.D.: 8.375 x 8.375 x 1.375

I.D.: 3.750 x 3.750 x 1.250

Mat.: OPS, HIPS



80065

4 Cell Tray

O.D.: 8.375 x 8.375 x 1.000

I.D.: 3.250 x 3.250 x 1.000

Mat.: OPS, HIPS



2401

4 Cell Tray

O.D.: 8.750 x 3.500 x 1.750

I.D.: 1.875 x 3.000 x 1.750

Mat.: OPS, HIPS



2440

4 Cell Tray

O.D.: 9.000 x 4.500 x 1.250

I.D.: 4.000 x 2.125 x 1.250

Mat.: OPS, HIPS



2400

4 Cell Tray

O.D.: 9.500 x 3.000 x 1.250

I.D.: 1.000 x 2.625 x 1.250

Mat.: OPS, HIPS



87204

4 Cell Tray

O.D.: 9.875 x 9.875 x 1.000

I.D.: 4.125 x 4.125 x 1.000

Mat.: PETE



**2405**

**5 Cell Tray**

**O.D.:** 11.125 x 3.750 x 1.250

**I.D.:** 1.500 x 3.000 x 1.250

**Mat.:** OPS, HIPS



**89966**

**6 cell Tray**

**O.D.:** 4.000 x 2.875 x 1.000

**I.D.:** 1.250 x 1.000 x 1.000

**Mat.:** OPS, HIPS



**88974**

**6 Cell Tray, Sealable**

**O.D.:** 4.125 x 3.000 x 0.500

**I.D.:** 1.125 x 1.125 x 0.500

**Mat.:** OPS, HIPS



**86506**

**6 Cell Tray**

**O.D.:** 7.125 x 4.250 x 1.250

**I.D.:** 1.875 x 1.875 x 1.250

**Mat.:** OPS, HIPS



**2508**

**6 Cell Tray**

**O.D.:** 8.000 x 5.375 x 1.750

**I.D.:** 2.125 x 2.125 x 1.750

**Mat.:** OPS, HIPS



**84906**

**6 Cell Tray**

**O.D.:** 10.500 x 4.750 x 1.250

**I.D.:** 3.250 x 2.000 x 1.250

**Mat.:** OPS, HIPS



**2990**

**6 Cell Tray**

**O.D.:** 10.625 x 7.250 x 1.625

**I.D.:** 3.125 x 3.125 x 1.625

**Mat.:** OPS, HIPS



**87421**

**6 Cell Tray**

**O.D.:** 13.000 x 8.500 x 1.000

**I.D.:** 3.375 x 3.375 x 1.000

**Mat.:** PETE



**89606**

**6 Cell Reverse Snap Tray**

**O.D.:** 13.250 x 8.750 x 2.000

**I.D.:** 3.250 x 3.250 x 0.375

**Mat.:** PETE



2807

6 Cell Tray, Stackable

O.D.: 19.750 x 13.125 x 2.375

I.D.: 5.875 x 5.875 x 2.375

Mat.: OPS, HIPS



88258

8 Cell Tray

O.D.: 7.625 x 5.125 x 0.750

I.D.: 1.875 x 1.875 x 0.750

Mat.: OPS, HIPS



2501

8 Cell Tray

O.D.: 9.750 x 5.000 x 1.375

I.D.: 2.125 x 2.125 x 1.375

Mat.: OPS, HIPS



2500

8 Cell Tray

O.D.: 10.250 x 5.250 x 1.240

I.D.: 2.250 x 2.250 x 1.000

Mat.: OPS, HIPS



2808

8 Cell Tray, Stackable

O.D.: 18.375 x 9.375 x 1.750

I.D.: 4.000 x 4.000 x 1.750

Mat.: OPS, HIPS



87424

10 Cell Tray

O.D.: 9.000 x 8.500 x 1.375

I.D.: 2.000 x 2.000 x 1.375

Mat.: OPS, HIPS



2120

12 Cell Tray

O.D.: 7.750 x 6.250 x 0.500

I.D.: 1.675 x 1.675 x 0.500

Mat.: OPS, HIPS




87228

12 Cell Tray

O.D.: 8.000 x 6.000 x 1.000

I.D.: 1.500 x 1.500 x 1.000

Mat.: OPS, HIPS



2123

12 Cell Tray

O.D.: 9.000 x 6.000 x 1.250

I.D.: 1.750 x 1.750 x 1.250

Mat.: OPS, HIPS





**81250**

**12 Cell Tray**

**O.D.:** 14.875 x 11.250 x 2.000

**I.D.:** 3.250 x 3.250 x 2.000

**Mat.:** OPS, HIPS



**81312**

**12 Cell Tray**

**O.D.:** 15.750 x 11.750 x 1.625

**I.D.:** 3.000 x 3.000 x 1.625

**Mat.:** OPS, HIPS



**87244**

**20 Cell Tray**

**O.D.:** 9.750 x 7.250 x 1.000

**I.D.:** 1.500 x 1.500 x 0.875

**Mat.:** PETE




**85205**

**15 Cell Tray**

**O.D.:** 9.375 x 6.000 x 0.875

**I.D.:** 2.000 x 2.000 x 0.875

**Mat.:** OPS, HIPS



**87215**

**15 Cell Tray**

**O.D.:** 9.563 x 6.000 x 0.875

**I.D.:** 1.375 x 1.375 x 0.875

**Mat.:** OPS, HIPS



**80619**

**16 Cell Tray**

**O.D.:** 8.625 x 8.625 x 1.188

**I.D.:** 2.000 x 2.000 x 1.125

**Mat.:** OPS, HIPS



**85196**

**16 Cell Tray**

**O.D.:** 8.625 x 8.625 x 0.750

**I.D.:** 1.875 x 1.875 x 0.750

**Mat.:** OPS, HIPS



**84820**

**18 Cell Tray**

**O.D.:** 10.625 x 6.000 x 0.875

**I.D.:** 1.125 x 1.125 x 0.625

**Mat.:** OPS, HIPS




**2040**

**1 Cell Mini Cheesecake Base**

**O.D.:** 3.750 x 3.750 x 0.500

**I.D.:** 2.750 x 2.750 x 0.250

**Mat.:** PETE



**2041**

**1 Cell Mini Cheesecake Lid**

**O.D.:** 3.750 x 3.750 x 2.044

**I.D.:** 3.500 x 3.500 x 1.675

**Mat.:** OPS, HIPS



**87781**

**Dome Lid, Fits PT 102**

**O.D.:** 4.375 x 4.375 x 2.125

**I.D.:** 3.125 x 3.125 x 2.000

**Mat.:** PETE



**5103**

**Dome Lid**

**O.D.:** 5.250 x 5.250 x 1.000

**I.D.:** 4.125 x 4.125 x 1.000

**Mat.:** OPS, HIPS



**2012**

**Lid, Snap On, Fits 6" Foam Bowl**

**O.D.:** 6.125 x 6.125 x 1.000

**I.D.:** 6.125 x 6.125 x 1.000

**Mat.:** OPS, HIPS



**88700**

**Dome Lid, Fits PT103**

**O.D.:** 8.000 x 5.375 x 1.500

**I.D.:** 7.500 x 4.125 x 1.500

**Mat.:** OPS, HIPS



**2018**

**Dome Lid**

**O.D.:** 8.500 x 6.500 x 0.875

**I.D.:** 7.500 x 5.125 x 0.875

**Mat.:** OPS, HIPS



**85031**

**1 Cell Base**

**O.D.:** 9.125 x 7.000 x 2.250

**I.D.:** 8.000 x 5.875 x 2.250

**Mat.:** PET




**85032**

**1 Cell Lid, Fits #85031**

**O.D.:** 9.375 x 7.250 x 0.500

**I.D.:** 9.125 x 7.000 x 0.500

**Mat.:** PET



**2013**

**Dome Lid, Fits PT109**

**O.D.:** 9.250 x 7.250 x 0.675

**I.D.:** 8.875 x 6.750 x 0.675

**Mat.:** OPS, HIPS



**2612**

**Dome Lid, Fits Waxboard Base**

**O.D.:** 11.000 x 6.000 x 3.125

**I.D.:** 10.000 x 4.500 x 3.000

**Mat.:** OPS, HIPS



**2019**

**Dome Lid, Crimp Over**

**O.D.:** 12.375 x 6.250 x 1.250

**I.D.:** 11.625 x 5.500 x 1.250

**Mat.:** OPS, HIPS




**87175**

**1/3 Tray Lid**

**O.D.:** 12.875 x 6.875 x 0.500

**I.D.:** 12.500 x 6.500 x 0.500

**Mat.:** PETE




**87502**

**1/4 Sheet Bakery Cover**

**O.D.:** 13.250 x 9.125 x 1.000

**I.D.:** 12.375 x 8.375 x 0.875

**Mat.:** OPS, HIPS



**87501**

**1/2 Sheet Bakery Cover**

**O.D.:** 18.250 x 13.375 x 0.750

**I.D.:** 17.250 x 12.375 x 0.625

**Mat.:** OPS, HIPS



**87246**

**Lid for 2 Cell Tray, Fits #87245**

**O.D.:** 10.250 x 5.250 x 0.625

**I.D.:** 10.000 x 5.000 x 0.500

**Mat.:** OPS, HIPS



**87245**

**2 Cell Tray**

**O.D.:** 10.000 x 5.000 x 1.000

**I.D.:** 4.125 x 4.125 x 1.000

**Mat.:** OPS, HIPS



**3009**

**1 Cell Hinge**

**O.D.:** 4.000 x 4.750 x 3.500

**I.D.:** 2.125 x 2.125 x 1.750

**Mat.:** PETE



**85024**

**1 Cell Hinge**

**O.D.:** 3.625 x 3.625 x 3.000

**I.D.:** 3.000 x 3.000 x 1.500

**Mat.:** PETE



87308

1 Cell Hinge

O.D.: 4.500 x 4.500 x 2.125

I.D.: 4.000 x 4.000 x 1.625

Mat.: PETE




2610

1 Cell Hinge w/Angled Back Edge

O.D.: 7.875 x 7.375 x 2.750

I.D.: 6.250 x 4.625 x 1.125

Mat.: PETE




86001

1 Cell Hinge

O.D.: 8.000 x 5.125 x 1.750

I.D.: 6.000 x 3.500 x 1.375

Mat.: PETE



2620

1 Cell Perimeter Seal Hinge

O.D.: 8.000 x 6.625 x 3.250

I.D.: 6.000 x 4.500 x 2.875

Mat.: PETE



2614

1 Cell Hinge

O.D.: 8.625 x 7.000 x 2.125

I.D.: 7.875 x 5.625 x 1.250

Mat.: PETE



80608

1 Cell Half Pie Hinge

O.D.: 8.750 x 5.250 x 1.750

I.D.: 7.675 x 3.875 x 1.125

Mat.: PETE



3434

1 Cell Hinge

O.D.: 8.875 x 7.875 x 2.750

I.D.: 7.000 x 6.000 x 1.250

Mat.: OPS, HIPS



87901

1 Cell Hinge

O.D.: 16.250 x 10.250 x 4.000

I.D.: 15.375 x 9.000 x 2.000

Mat.: PETE



3010

4 Cell Hinge

I.D.: 4.750 x 4.750 x 2.125

O.D.: 1.625 x 1.625 x 1.125

Mat.: PETE





**2613**

**4 Cell Hinge**

**O.D.:** 8.000 x 7.500 x 3.375

**I.D.:** 2.625 x 2.625 x 1.625

**Mat.:** PETE



**89802**

**4 Cell Hinge**

**O.D.:** 8.000 x 7.500 x 3.000

**I.D.:** 2.625 x 2.625 x 1.625

**Mat.:** PETE



**86806**

**6 Cell Hinge**

**O.D.:** 7.125 x 4.750 x 1.875

**I.D.:** 1.375 x 1.375 x 1.750

**Mat.:** PETE



**2615**

**6 Cell Hinge**

**O.D.:** 8.625 x 7.000 x 2.125

**I.D.:** 2.125 x 2.125 x 1.250

**Mat.:** PETE



**84358**

**8 Cell Hinge**

**O.D.:** 9.250 x 4.750 x 2.000

**I.D.:** 1.875 x 1.875 x 1.000

**Mat.:** OPS, HIPS



**84359**

**10 Cell Hinge**

**O.D.:** 11.625 x 6.125 x 2.250

**I.D.:** 1.750 x 2.500 x 1.125

**Mat.:** OPS, HIPS



**2611**

**12 Cell Hinge**

**O.D.:** 8.675 x 7.000 x 2.125

**I.D.:** 2.000 x 2.000 x 1.250

**Mat.:** PETE



**85027**

**12 Cell Hinge**

**O.D.:** 10.000 x 12.250 x 2.750

**I.D.:** 2.125 x 2.125 x 1.500

**Mat.:** PETE